

# Guidelines for Keratoconus Data Collection

## Instructions to the Patient

1. Ask the patient to **close** one-eye with their **hand**.
2. Ask the patient to **blink a few times** and **keep the other eye wide open**.
3. Ask the patient to **look at the camera (center of the rings)**.

## Instructions to the Operator

During data collection, please ensure:

1. Correctly enter the **patient id**, **age** and **gender**.
2. Ensure that the **LED lights are ON** (blue color will appear on the screen).
3. Place the attachment around the patient's eye orbit.
4. Both the **top and bottom part of the attachment should touch the face**.
5. The ring reflections should be **centered**.
6. The camera is **in focus**. Please tap to focus, if not.
7. Reflection image (mires) should cover **limbus to limbus**.
  - a. NO minified images.
8. The two side **supports should be horizontal**.

## To select an image (out of the three clicked images):

1. The image is **in focus, with no blur**.
2. The two side **supports are horizontal**.
3. The ring reflections are **centered**.
4. Good **limbus to limbus** coverage.

### Note:

- At the end of each day, upload the reports (mire image, axial image, and/or tangential image) from Keratron and/or Occulyzer device.
- Note: Internet connection is required for SmartKC App to upload data.

## For each patient, perform the following steps:


1. Optikon Keratron or Occulyzer Test.
2. Collect data from the phone app with placido attachment.

## Good and Bad Image Examples

If the above guidelines are not followed properly, it may result in a bad image capture that can lead to incorrectly generated topography maps and a misleading diagnosis. Hence, it is of utmost importance to follow guidelines properly at all times. Below are examples of good and bad image scenarios, that can act as a reference when training an operator to use the SmartKC device to capture perfect images:



### Good Image Examples

1. Mires sharp and in-focus, limbus-to-limbus coverage, support horizontal, no offset





### Bad Image Examples

1. Bad limbus-to-limbus coverage (minified image), outside reflections visible (placido head not touching area around eye orbital well)



External Reflections

2. Mires blurred, bad limbus-to-limbus coverage, external reflections visible.





Mires blurred

3. Mires blurred due to **tear-film** breaking, and also has bad limbus-to-limbus coverage.



4. Mires broken due to eye-lashes

Eyelashes break mires





5. Mires bleeding out of iris region into the sclera.

